

WARNING NOTICE: The experiments described in these materials are potentially hazardous and require a high level of safety training, special facilities and equipment, and supervision by appropriate individuals. You bear the sole responsibility, liability, and risk for the implementation of such safety procedures and measures. MIT shall have no responsibility, liability, or risk for the content or implementation of any of the material presented. [Legal Notices](#)

ELECTROCOMPETENT *RHODOCOCCUS* I24 AND KY1

1. (Day or two before) Inoculate small (5-10 ml) cultures of *Rhodococcus* I24; grow at 37°C
2. Transfer 0.1-5 ml of overnight culture of *Rhodococcus* to 200 ml MB 1.5% Glycine (for I24 strain) or TSB Glycine (for KY1 strain) in a 1L baffled flask
3. Incubate shaking at 37°C overnight or until O.D.₆₀₀ is approximately 0.25
4. Pellet the cells by centrifuging for 5 min at 6 000 rpm in a GSA rotor using sterile centrifuge bottles or 50 ml conical tubes and proper adapters (may have to spin twice to pool)
5. Resuspend the cell pellet in 30 ml ice-cold EPB1; Recentrifuge as in step 4
6. Wash cell pellet one more time in EPB1; centrifuge as before; discard supernatant
7. Wash pellet once in 10 ml ice-cold EPB2; centrifuge as before except at 8000 rpm; discard supernatant
8. Resuspend final cell pellet in 1 ml or less of EPB2
9. Aliquot 150 µl into sterile microfuge tubes and store at -80°C

Electroporation of *Rhodococcus*

7. Thaw aliquots of electrocompetent *Rhodococcus* cells on ice
8. Mix DNA with 70µl cells in a sterile microfuge tube and incubate on ice for 5 min.
9. Electroporate DNA at 2.5 kV, 25 µF and 400 Ω
10. Immediately add 300 µl LB
11. Incubate cells for recovery at 30°C for 1-20 hours
12. Spread cells onto plates with appropriate antibiotics

MB 1.5% Glycine medium (per liter)

Yeast extract	5g
Bacto tryptone	15 g
Bacto soytone	5g
NaCl	5g
Glycine	15g

TSB Glycine(per liter)

Per liter	
Bacto Tryptone	17 g
Bacto Soytone	3 g
Sodium Chloride	5 g
Dipotassium Phosphate	2.5 g
Yeast Extract	5 g
Glycine	15 g

Hepes Stock Solution

Hepes	23.8g
distilled water	180ml
adjust pH to 7.2; raise volume to 200 ml	

EPB1 (20 mM Hepes, 5% glycerol, pH7.2)

0.5 M Hepes stock, pH7.2	20ml
100% glycerol	25ml
distilled water to 500 ml	

EPB2 (5mM Hepes, 15% glycerol, pH7.2)

0.5 M Hepes stock, pH7.2	2ml
100% glycerol	30ml
distilled water to 200ml	