

WARNING NOTICE: The experiments described in these materials are potentially hazardous and require a high level of safety training, special facilities and equipment, and supervision by appropriate individuals. You bear the sole responsibility, liability, and risk for the implementation of such safety procedures and measures. MIT shall have no responsibility, liability, or risk for the content or implementation of any of the material presented. [Legal Notices](#)

QUANTITATING NUCLEIC ACIDS

DNA or RNA (based on absorbance at 260 nm)

1. prepare a dilution of DNA or RNA sample in distilled water (typically 1:100 or 5:500 μ l)
2. using UV light source, read absorbance at 260nm, may also check absorbance at 280 and 230nm (260/280 ratio should be roughly 2/1; 260/230 ratio should also approach 2/1)
3. to calculate concentration of original solution:

$A_{260} \times \text{conversion factor} \times \text{dilution factor} = \text{original DNA concentration } (\mu\text{g/ml})$ absorbance unit
the conversion factor for double stranded DNA is 50 $\mu\text{g/ml}$; for single stranded DNA and RNA
the conversion factor is 40 $\mu\text{g/ml}$

The molecular weight of an average deoxynucleotide monophosphate (base) is 326.95. Therefore, the average basepair is 753.9. To calculate μ moles of ends from this value, for example, calculate the following:

$$(2 \times \mu\text{g DNA}) + [753.9 \times (\text{size of double stranded DNA in bp})]$$

thus 1 μg of pUC DNA corresponds to $2 \times 1\mu\text{g} + [753.9 \mu\text{g}/\mu\text{mol bp} \times 2676 \text{ bp}]$ or $2\mu\text{g} + 2017436.4 \mu\text{g}/\mu\text{mol} = 1 \times 10^{-6}\mu\text{mol}$

Oligodeoxynucleotides

1. calculate extinction coefficient (ϵ) for oligonucleotide:
 - a. tabulate number of occurrences of each base, A, C, G and T
 - b. multiply results by the following values:

number of A's x	15.4 ml/ μ mole
number of C's x	7.3 ml/ μ mole
number of G's x	11.7 ml/ μ mole
number of T's x	<u>8.8 ml/μmole</u>
 - c. the sum of these results is the extinction coefficient, ϵ
2. determine absorbance of oligonucleotide solution at 260 nm
3. calculate concentration with the following formula:

$$\text{concentration} = A_{260} \div \epsilon \text{ (note: } \mu\text{mole/ml} = \text{mM)}$$

NB. oligonucleotides are frequent supplied as lyophilized powders with amounts listed in O.D. units; one O.D. unit in this nomenclature corresponds to the A_{260} this sample would have were it resuspended in 1ml (roughly); this information can be used with the above formulae to estimate what the concentration of oligo would be in such a solution and a more preferred concentration of oligo can be prepared based on this information;

for example:

an oligo has a ϵ of 262.9 ml/ μ mole; 3 O.D. units were supplied; therefore, a 1ml solution would have a concentration of $3 \div 262.9 \text{ ml}/\mu\text{mole} = 11.4 \mu\text{M}$ to obtain a, say, 20 μM solution calculate $C_1V_1 = C_2V_2$:

$$1\text{ml} \times 11.4 \mu\text{M} = ? \times 20 \mu\text{M}; \quad (1\text{ml} \times 11.4\mu\text{M})/20\mu\text{M} = 0.570\text{ml}$$

resuspend in 570 μ l!