

WARNING NOTICE: The experiments described in these materials are potentially hazardous and require a high level of safety training, special facilities and equipment, and supervision by appropriate individuals. You bear the sole responsibility, liability, and risk for the implementation of such safety procedures and measures. MIT shall have no responsibility, liability, or risk for the content or implementation of any of the material presented. [Legal Notices](#)

SOUTHERN BLOT PROTOCOL USING DIGOXIGENIN LABELED PROBE, RHODOCOCCUS OPTIMIZED

selflessly tested , optimized and illustrated with love by xian o'brien

Recipes for reagents in **bold** included at the end of this protocol, in *italics* are included with the Roche DIG labeling and detection kit

LABELING DNA PROBE USING DIG HIGH PRIME LABELING MIX (ROCHE)

- dilute 10ng-3 μ g probe DNA (genomic, plasmid or gene clean fragment) in dsH₂O to a final volume of 16 μ l
- denature DNA by boiling 10min; quickly chill on ice to prevent reannealing of strands
- add 4 μ l *DIG high prime labeling mix*; mix briefly and tap spin
- incubate overnight at 37°C.
- stop reaction by adding 2 μ l 0.2M EDTA (pH 8) and heat inactivate at 65°C 10min
- boil probe 10-20min before using
- used probe can be stored at -20°C in **Hybridization Buffer** and used repeatedly

TRANSFER OF DNA FROM GEL TO MEMBRANE


- run gel at appropriate voltage to obtain clean bands; stain and photograph with size reference
- smaller fragments transfer more efficiently. for very large DNAs, a 2min UV nicking step on a short wave transilluminator can be added
- transfer gel to a sealable Tupperware container

traditional method of transfer

- incubate 40min at room temperature in **0.25 HCl** (sufficient to cover gel) to depurinate
- rinse 2X in MilliQ
- incubate 2X 20min in **Denaturation Solution** to cleave depurination sites
- incubate 30min in **Neutralization Solution**
- set up capillary transfer as shown in the schematic below using **20X SSC** as the Transfer Solution
- transfer overnight (48hr for pulsed field gels)

alternative quick transfer method (Phil's favorite)

- incubate at room temperature in **0.25M HCl** until dye bands turn yellow (*ca.* 20min)
- rinse 2X in MilliQ
- set up capillary transfer as shown in the schematic below using **0.4N NaOH** as the Transfer Solution
- transfer overnight (48hr for pulsed field gels)



SOUTHERN BLOT (CONT.)

FIX DNA TO MEMBRANE

- cut off left top corner of gel and membrane for orienting blot
- restain and photograph gel to assess transfer efficiency
- rinse membrane briefly in **2X SSC** and transfer to sealable Tupperware container
- optional*: prestrip membrane using **BLOT STRIPPING PROTOCOL**. often yields cleaner blots

PROBE MEMBRANE

- prehyb membrane on rocker for a minimum of 3hr at 68°C in **Hybridization Buffer**
- boil probe in 40ml **Hybridization Buffer** at least 10min to denature
- hybridize membrane DNA side down overnight at 68°C (42° for lower homology) in boiled **Hybridization Buffer/Probe mix**, rocking optional
- used **Hybridization Buffer/Probe mix** can be stored at -20°C and used repeatedly

MEMBRANE DETECTION

- wash 2X 15min at room temperature with **2X SSC, 0.1% SDS** on rocker
- wash 2X 15min at 42°C with **0.5X SSC, 0.1% SDS** on rocker (washes can be modified to control stringency – this is fairly stringent)
- rinse in **Washing Buffer**
- incubate 30min at room temperature in **Blocking Solution**, rocking optional
- incubate 30min at room temperature with 30ml **Blocking Solution** + 2µl *anti-DIG-AP Conjugate* (premix before adding to blot), rocking optional
- wash 2X 15min at room temperature with 100ml **Washing Buffer** on rocker
- equilibrate 5min at room temperature with **Detection Buffer**
- lay membrane on saran wrap; add 20 drops *Ready-To-Use CSPD Reagent*
- cover with a second piece of saran wrap
- let stand 3min; squeeze out excess *CSPD Reagent* from between sheets of plastic wrap remove as much as possible to ensure low background on film
- expose to film (enzymatic reaction is accelerated at 37°C and slowed at 4°C)
- after developing, membranes can be stored as is between sheets of plastic wrap at room temperature indefinitely

MEMBRANE STRIPPING

- remove membrane from plastic wrap and place in a sealable Tupperware container
- wash 2X 15min at 37°C in **Stripping Buffer** (no longer!)
- rinse briefly in **2X SSC**
- membrane is now ready to prehyb

SOUTHERN BLOT (CONT.)

REAGENTS

0.25N HCl
25ml conc. HCl
MilliQ to 1L

Denaturation Solution
88g NaCl (1.5M)
20g NaOH (0.5N)
MilliQ to 1L

Neutralization Solution
176g NaCl (3M)
6.7g Tris base (0.5M)
70.2g Tris·HCl
MilliQ to 1L

20X SSC
176g NaCl (3M)
88g Na₃Citrate (0.3M)
MilliQ to 1L
pH to 7.0 w/ 1M HCl

2X SSC
50ml 20X SSC
MilliQ to 500ml

10X Maleic acid Buffer
116g maleic acid (1X is 0.1M)
88g NaCl (1X is 0.15M)
MilliQ to 1L
PH to 7.5 w/ solid NaOH

2X SSC; 0.1% SDS
100ml 20X SSC
10ml 10% SDS
MilliQ to 1L

0.5X SSC; 0.1% SDS
25ml 20X SSC
10ml 10% SDS
MilliQ to 1L

10X Detection Buffer
1M Tris·HCl
1M NaCl
pH to 9.5

Washing Buffer
3ml Tween-20 (0.3%)
100ml 10X maleic acid buffer
MilliQ to 1L

Blocking Solution
50ml 10X Blocking Reagent*
50ml 10X maleic acid buffer
MilliQ to 1L

Hybridization Buffer
250ml 20X SSC (5X)
100ml 1% lauryl sarcosine (0.1%)
2ml 10% SDS (0.02%)
100ml 10X Blocking Reagent*
MilliQ to 1L

10N NaOH
200g NaOH
MilliQ to 500ml

0.4N NaOH
40ml 10N NaOH
MilliQ to 1L

Stripping Buffer
8g NaOH (0.2N)
10ml 10% SDS (0.1%)
MilliQ to 1L

* included in Roche DIG kit