

Course 12.141: Electron Microprobe Analysis – Problem Set 2

Problem 2.1: Use wavelength dispersive spectrometry (WDS) to acquire X-ray spectra to identify elements in each phase of the sample provided. Perform semi-quantitative analysis to obtain approximate composition of the phases.

Problem 2.2: Acquire X-ray elemental maps of the polished sample provided with WDS and EDS. Identify the phases and use an image analysis software to find the area percent of each phase.

MIT OpenCourseWare
<http://ocw.mit.edu>

12.141 Electron Microprobe Analysis
January (IAP) 2012

For information about citing these materials or our Terms of Use, visit: <http://ocw.mit.edu/terms>.